

# ChIP protocol – preparation from liver

Version 3.3 – CIGENE/NMBU

## Day 1 – Tissue disruption, crosslinking and sonication

### Preparations

- Equipment:
  - microscope, cell counting chamber, Trypan Blue/Hoescht
  - tweezers
  - a dounce homogenizer with pestles <sup>1</sup>
  - a probe sonicator / sonication device
  - an ice bucket
  - a 15-ml centrifuge set-up at 4°C
  - a 1.5/2-ml centrifuge set-up at 4°C
  - a metallic eppendorf tubes holder
  - an incubator set-up at 68°C
- Solutions (see **Buffers for ChIP**):
  - PBS with PIC
  - formaldehyde solution 1% (keep in fridge until use)
  - sonication buffer completed with PIC, PMSF and sodium butyrate (named **complete sonication buffer** in the protocol)
  - IP buffer completed with PIC, PMSF and sodium butyrate (named **complete IP buffer** in the protocol)
- To do:
  - refrigerate the metallic Eppendorf tubes holder in the fridge at 4°C
  - Pour 4.5 ml of **PBS with PIC** in the dounce homogenizer, and place it in the ice bucket
  - take the appropriate volume of **formaldehyde solution 1%** from the fridge and place it at room temperature before use

### Tissue disruption

**The procedure (step 1 to 3) should take less than 15 min and the samples should stay in cold environment as much as possible. Glass homogenizers must be used in priority to plastic ones, due to better conductivity of cold**

- 1** Transfer pieces of liver (60 mg recommended for liver) into the dounce homogenizer containing 4.5 ml <sup>2</sup> of **PBS with PIC**
- 2** Use pestles A and B for the homogenization, 5 strokes each or until obtaining a homogenized solution or feel no more strong resistance
- 3** Transfer the homogenate into a 5-ml eppendorf tube <sup>3</sup>. You can cut the bottom of the pipet tip to facilitate the aspiration
- 4** Rinse the douncer and the pestles with 0.5 ml of **PBS with PIC**. Add it to the previous solution. You reach 5ml
- 5** After homogenization of the tube, take an aliquot of **10 µl to assess the number of nuclei** (dilute 10x, add trypan blue or Hoescht for microscopy observation)  
**Counting:**                      **nuclei/ml**                      **5ml =**                      **million nuclei total** (expected: 20 million)
- 6** Centrifuge 3 min, 2 500 g, 4°C
- 7** Remove and discard supernatant

### Crosslinking

- 8** Add to the pellet 2.5 ml\* of **formaldehyde solution 1%** at room temperature (RT) and resuspend the pellet by pipetting  
*\*if more than 24 million of nuclei, resuspend in 5 ml of formaldehyde solution 1% and split in two tubes*
- 9** Incubate under constant agitation <sup>4</sup> at room temperature, 5 min
- 10** Quench the reaction with 360 µl of **glycine 1M** (0.125M final). Incubate 10 min at RT, under constant agitation <sup>2</sup>
- 11** Centrifuge 5 min, 800 g, 4°C
- 12** Remove and eliminate in appropriate trash the supernatant containing formaldehyde
- 13** Resuspend with 4 ml\* of **PBS with PIC** by pipetting (quick vortex if necessary)

<sup>1</sup> 2-ml or 15-ml dounce

<sup>2</sup> for a 2-ml dounce, homogenise in 1.5ml solution. Then rinse the dounce with the rest of the solution and transfer into the 5-ml Eppendorf tube

<sup>3</sup> or a 15-ml canonical Falcon tube

<sup>4</sup> Place the tube in the Hula Mixer, 360° rotation, 40 rpm. Other options : in horizontal position in a box on the top of an orbital shaker/rocker (60 rpm) or directly between two tubes of a tube roller mixer (60 rpm)



## Version history

<b>Version 3.3 – 01.10.2020</b>	Comment about sonication profile change
<b>Version 3.2 – 01.10.2020</b>	Add of a version history section Change of the safe stop point explanation