

ChIP protocol – all tissues

Version 3.2 – CIGENE/NMBU

Day 3 – Wash beads, decrosslinking and DNA recovery

Preparations

- Equipment:
 - a 1.5/2ml centrifuge set-up at 4°C
 - a magnetic rack
 - a mini-centrifuge
 - an incubator set-up at 65°C/68°C
- Solutions (see **Buffers for ChIP**):
 - IP buffer completed with PIC, PMSF and sodium butyrate (named **complete IP buffer** in the protocol)
 - Low salt wash buffer completed with PIC, PMSF and sodium butyrate (named **low salt wash buffer** in the protocol)
 - High salt wash buffer completed with PIC, PMSF and sodium butyrate (named **high salt wash buffer** in the protocol)
 - ChIP elution buffer and elution buffer

Wash of the beads

- 1 Pulse spin the tubes, put on magnet and discard supernatant
 - 2 Resuspend the beads in 200 µl of **complete IP buffer**
 - 3 Pulse spin the tubes, put on magnet and discard supernatant
 - 4 Resuspend the beads in 200 µl of **low salt wash buffer**, put in the magnet, discard supernatant
 - 5 Repeat **step 4**
 - 6 Resuspend the beads in 200 µl of **high salt wash buffer**, put in the magnet, discard supernatant
 - 7 Repeat **step 6**
 - 8 Resuspend the beads in 70 µl of **ChIP elution buffer**
 - 9 Elute the DNA from the beads for 1H30 at 65°C, under agitation ¹
 - 10 Pulse spin the tubes and place on magnetic rack
 - 11 Transfer the eluted chromatin to new clean LoBind 1.5-ml Eppendorf tubes
 - 12 Wash the beads with 30 µl of **ChIP elution buffer**. Mix and pulse spin the tubes
 - 13 Place the tubes on magnetic rack. Transfer this volume to the eluted chromatin from **step 11** (100 µl final)
 - 14 Take out the "input control" tube stocked in the fridge
 - 15 Complete the volume of the input to 100 µl final
 - 16 Add 5 µl of SDS 20% (1% final of SDS) into the input tube, mix well
- The input tube and the eluted chromatin from IP tubes are ready for DNA extraction**

One-step decrosslinking

- 17 To each tube, add 2 µl of **RNase A** (100 mg/ml), 5 µl of **proteinase K** (20 mg/ml) and 6 µl of **NaCl 5M**
- 18 Incubate 1H30, 68°C, under agitation ²
- 19 Take out from the fridge the **Qubit HS kit**

Purification and quantification

- 20 Purify the samples with **MinElute PCR purification kit**. Elute from the column with 32 µl of **elution buffer**
- 21 Quantify 2 µl of the samples with **Qubit HS kit**

The DNA from the input and immunoprecipitated chromatin can be stored at – 20°C until qPCR and library preparation

¹ 1000 rpm into an Eppendorf Thermomixer

² 500 rpm into an Eppendorf Thermomixer

Version 3.2 – 01.10.2020	Add of a version history section
---------------------------------	----------------------------------